

Live Cell Imaging Resource Lab




Practical approaches to optical microscopy education at the LCI resource lab

RIMA WAZEN, PHD

LIVE CELL IMAGING RESOURCE LAB

UNIVERSITY OF CALGARY


Knowledge gap in microscopy



What are the options for hands-on microscopy training?

- ❖ Montreal Light Microscopy Course (MLMC), Montreal, QC
- ❖ Canadian Light Microscopy Course (CLMC), Calgary, AB
- ❖ Marine Biological Laboratory, Woods Hole, MA
- ❖ Cold Spring Harbor Laboratory, Cold Spring Harbor, NY
- ❖ Mount Desert Island (MDI) Biological Laboratory, Bar Harbor, ME
- ❖ Institut Pasteur, Paris, France
- ❖ European Molecular Biology Laboratory (EMBL), Heidelberg, Germany
- ❖ Workshops offered by commercial organizations
- ❖ and many more...

How are we helping our research community?



Mission: Drive research, education and training in optical imaging

Educational activities at all level

- ❖ Step-wise training program for the LCI's users
- ❖ Workshops and educational opportunities for researchers who do not use the LCI instruments

Step-wise training program to the LCI users

- ❖ Meeting with the researcher to **define the project** imaging goal
 - Experimental design
 - Image acquisition
 - Image analysis plan
- ❖ Pre-training **worksheet** and review **learning outcomes**
 - Introduction to fluorescence & wide-field microscopy
 - Laser scanning confocal microscopy
 - Spinning disk confocal microscopy
 - Multiphoton microscopy
 - Image analysis
- ❖ One-on-one training session(s) tailored to the project and the user's background
- ❖ Follow up meetings and on-going support with the researchers




Consolidating the skill set



Monthly (or as needed) microscopy workshops

- *Brightfield imaging*
- *White-light and contrast techniques*
- *Fluorescence wide-field microscopy*



Sample Preparation & Immunofluorescence



ImageJ
Image Processing & Analysis in Java

Image analysis

- *Digital imaging analysis*
- *Colour and colocalization*
- *Ethics of digital imaging*

Educational opportunities



UNIVERSITY OF
CALGARY

*3 Graduate
Microscopy Modules*



*3-day workshop
Fundamentals of microscopy*



Canadian Light Microscopy Course
July 17-21, 2017

*5-day workshop
Intermediate to advanced
microscopy*

Example of the curriculum for the 3-day LCI jumpstart workshop



*Microscope on a rail:
Understanding image formation*



*White-light &
Contrasting techniques*





*Intro to fluorescence microscopy
& Digital imaging acquisition*

Over the years, the following people have contributed to the jumpstart:
*Craig Brideau, Claire Brown, Andrew Chojnacki, Pina Colarusso, Grant Gordon,
Jennifer Poirier, Kelvin Poon, Rima Wazen, Katarzyna Wojcik, and many others*

Example of the curriculum for the 3-day LCI jumpstart workshop

Example of application-based activity

- *Acquisition of multi-color images*
- *Spatial and temporal sampling*
- *Setting up for live cell imaging*



HeLa cells stably expressing MitoDsRed.
Courtesy of Nicole Mancini, McKay lab

Example of the curriculum for the 3-day LCI jumpstart workshop

➤ **Case studies:** bridging the gap between theory and good practice

John needs multi-channel images of some fixed samples. He wants to understand the spatial relationships between proteins A, B, and C at the highest resolution possible. Proteins A and B are abundant and label well, but the protein C labelling is dim and hard to see. He has two microscope systems available to him, one with a sCMOS camera (QE=80%) and the other with a CCD (QE=55%). Both microscopes are adjusted to image at Nyquist at 40X magnification as purchased, but the CCD's photosites are twice the size of the sCMOS sensors. He has a 40X 1.3 NA for the CCD system and a 40X 0.75NA objective for the sCMOS system. Both microscopes have identical filter sets (GFP/RFP/Cy5) and identical light sources. You are John's PI. How would you advise John to label his proteins and which scope should he use?



Example of the curriculum for the 3-day LCI jumpstart workshop

Best practices for instrument maintenance

- ❖ Cleaning objectives & filters, Check dust
- ❖ Inspecting even illumination in epifluorescence microscopy
- ❖ Checking system behavior (drift test, multi-point visiting, stitching, multi-color, etc...)


Future directions

- ❖ Create partnerships with optical imaging centres across Canada to strengthen, refine and distribute the LCI's educational and training programme.

- ❖ Resource package for imaging platform manager including
 - Education material for workshops and community outreach events
 - Pre-training worksheets
 - Best practices for instrument maintenance
 - Training records & user agreement

Acknowledgments

- Andrew Chojnacki, PhD, In vivo Optical Imaging Specialist
- Pina Colarusso, PhD, Director
- Joel Glover, Optical imaging assistant
- Kamala Patel, PhD, Executive Director
- Katarzyna Wojcik, PhD, In vivo Optical Imaging Specialist



UNIVERSITY OF CALGARY
CUMMING SCHOOL OF MEDICINE